

4E-BP1 / EIF4EBP1 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological
Biological Solution Specialist

Catalog Number: 10022-T16

GENERAL INFORMATION

Immunogen:	Recombinant Human 4E-BP1 / EIF4EBP1 protein (Catalog#10022-H07E)
Preparation	Produced in rabbits immunized with purified, recombinant Human 4E-BP1 / EIF4EBP1 (rh 4E-BP1 / EIF4EBP1; Catalog#10022-H07E; NP_004086.1; Ser2-Ile118). 4E-BP1 / EIF4EBP1 specific IgG was purified by Human 4E-BP1 / EIF4EBP1 affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Human 4E-BP1 / EIF4EBP1
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA
----------------------	-------

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human 4E-BP1 / EIF4EBP1.
--------------	-------------------------------------------------------------------------------------------------------------------------------------------------

Please Note: Optimal concentrations/dilutions should be determined by the end user.