

2B4 / SLAMF4 / CD244 Antibody, Mouse MAb



Sino Biological
Biological Solution Specialist

Catalog Number: 10042-MM02

GENERAL INFORMATION

Immunogen:	Recombinant Human 2B4 / SLAMF4 / CD244 Protein (Catalog#10042-H08H)
Preparation	This antibody was produced from a hybridoma resulting from the fusion of a mouse myeloma with B cells obtained from a mouse immunized with purified, recombinant Human 2B4 / SLAMF4 / CD244 (rh 2B4 / SLAMF4 / CD244; Catalog#10042-H08H; NP_057466.1; Met 1-Arg 221). The IgG fraction of the cell culture supernatant was purified by Protein A affinity chromatography.
Ig Type:	Mouse IgG1
Clone ID:	3F3A6
Specificity:	Human 2B4 / SLAMF4 / CD244
	No cross-reactivity in ELISA with Human SLAMF1 Human cell lysate (293 cell line)
Formulation:	0.2 µm filtered solution in PBS with 5% trehalose
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: ELISA

RECOMMENDED CONCENTRATION

ELISA ELISA: 1:1000-1:2000
This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Human 2B4.

Please Note: Optimal concentrations/dilutions should be determined by the end user.