



Catalog Number: 10042-R020

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION	
<b>Immunogen:</b>	Recombinant Human 2B4 / CD244 protein (Catalog#10042-H08H)
<b>Preparation</b>	This antibody was obtained from a rabbit immunized with purified, recombinant Human 2B4 / CD244 (rh 2B4 / CD244; Catalog#10042-H08H; NP_057466.1; Met 1-Arg 221).
<b>Ig Type:</b>	Rabbit IgG
<b>Clone ID:</b>	020
<b>Specificity:</b>	Human 2B4 / CD244
<b>Formulation:</b>	0.2 µm filtered solution in PBS with 5% trehalose
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
<b>Applications:</b>	ELISA
RECOMMENDED CONCENTRATION	
<b>ELISA</b>	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Human 2B4.

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**