



Catalog Number: 10042-R031

EliteRmab® is a registered trademark of Sino Biological Inc.

GENERAL INFORMATION

Immunogen:	Recombinant Human 2B4 / CD244 Protein (Catalog#10042-H08H)
Preparation	This antibody was obtained from a rabbit immunized with purified, recombinant Human 2B4 / CD244 (rh 2B4 / CD244; Catalog#10042-H08H; NP_057466.1; Met1-Arg221).
Ig Type:	Rabbit IgG
Clone ID:	031
Specificity:	Human 2B4 / CD244
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA(Det)
----------------------	------------

RECOMMENDED CONCENTRATION

Sandwich ELISA (Detection Ab)	ELISA(Det): 1:1000-1:10000 This antibody will detect Human 2B4 / CD244 in ELISA pair set (Catalog: # SEK10042). In a sandwich ELISA, it can be used as detection antibody when paired with (Catalog: # 10042-R025).
--------------------------------------	--

Please Note: Optimal concentrations/dilutions should be determined by the end user.