

**Influenza A H5N1 (A/Hong Kong/483/97)
Hemagglutinin / HA Antibody, Rabbit PAb,
Antigen Affinity Purified**

Catalog Number: 11689-T54



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION	
Immunogen:	Recombinant Influenza A H5N1 (A/Hong Kong/483/97) Hemagglutinin / HA Protein (Catalog#11689-V08H)
Preparation	Produced in rabbits immunized with purified, recombinant Influenza A H5N1 (A/Hong Kong/483/97) Hemagglutinin / HA (Catalog#11689-V08H; AAC32099.1; Met1-Gln531). Influenza A H5N1 (A/Hong Kong/483/97) Hemagglutinin / HA specific IgG was purified by Influenza A H5N1 (A/Hong Kong/483/97) Hemagglutinin / HA affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Influenza A H5N1 (A/Hong Kong/483/97) Hemagglutinin / HA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
Applications:	WB,ELISA,ICC/IF,IF,IP (Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)
RECOMMENDED CONCENTRATION	
Western Blot	WB: 1:1000-1:5000
ELISA	ELISA: 1:5000-1:10000 This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Influenza A H5N1 (A/Hong Kong/483/97) Hemagglutinin / HA.


Please Note: Optimal concentrations/dilutions should be determined by the end user.

**Influenza A H5N1 (A/Hong Kong/483/97)
Hemagglutinin / HA Antibody, Rabbit PAb,
Antigen Affinity Purified**

Catalog Number: 11689-T54



Sino Biological
Biological Solution Specialist



Anti-Influenza A H5N1 (A/Hong Kong/483/97)
Hemagglutinin / HA polyclonal antibody at
1:1000 dilution.

Sample: Influenza A H5N1 (A/Hong
Kong/483/97) Hemagglutinin / HA
Recombinant Protein

Lane A: 50ng

Lane B: 10ng

Secondary

Goat Anti- Rabbit IgG H&L (Dylight 800) at
1/10000 dilution.

Developed using the Odyssey technique.

Performed under reducing conditions.