

Influenza A H5N3 Hemagglutinin / HA Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11696-RP01

GENERAL INFORMATION

Immunogen:	Recombinant H5N3 (A/duck/Hokkaido/167/2007) HA / Hemagglutinin protein (Catalog#11696-V08H)
Preparation	Produced in rabbits immunized with purified, recombinant H5N3 (A/duck/Hokkaido/167/2007) HA / Hemagglutinin (Catalog#11696-V08H; BAG07130.2; Met 1-Gln 527). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	H5N3 (A/duck/Hokkaido/167/2007) HA / Hemagglutinin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications: WB,ELISA

RECOMMENDED CONCENTRATION

Western Blot	This antibody can be used at 1:500-1:1000 with the appropriate secondary reagents to detect H5N3 (A/duck/Hokkaido/167/2007) HA in WB.
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect H5N3 (A/duck/Hokkaido/167/2007) HA.

Please Note: Optimal concentrations/dilutions should be determined by the end user.