

Influenza A H13N8 Hemagglutinin / HA Antibody, Rabbit PAb



Sino Biological
Biological Solution Specialist

Catalog Number: 11721-RP01

GENERAL INFORMATION

Immunogen:	Recombinant H13N8 (A/black-headed gull/Netherlands/1/00) HA / Hemagglutinin protein (Catalog#11721-V08H)
Preparation	Produced in rabbits immunized with purified, recombinant H13N8 (A/black-headed gull/Netherlands/1/00) HA / Hemagglutinin (Catalog#11721-V08H; AAV91212.1; Met 1-Lys 528). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	H13N8 (A/black-headed gull/Netherlands/1/00) HA / Hemagglutinin
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	WB,ELISA
----------------------	----------

RECOMMENDED CONCENTRATION

Western Blot	This antibody can be used at 1:500-1:1000 with the appropriate secondary reagents to detect H13N8 (A/black-headed gull/Netherlands/1/00) HA in WB.
ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect H13N8 (A/black-headed gull/Netherlands/1/00) HA.

Please Note: Optimal concentrations/dilutions should be determined by the end user.