

Influenza A H7N9 (A/Shanghai/1/2013) Hemagglutinin / HA Antibody, Rabbit PAb

Catalog Number: 40104-RP01



Sino Biological
Biological Solution Specialist

GENERAL INFORMATION

Immunogen:	Recombinant Influenza A H7N9 (A/Shanghai/1/2013) Hemagglutinin / HA Protein (Catalog#40104-V08B)
Preparation	Produced in rabbits immunized with purified, recombinant Influenza A H7N9 (A/Shanghai/1/2013) Hemagglutinin / HA (Catalog#40104-V08B; Met1-Val524). Total IgG was purified by Protein A affinity chromatography.
Ig Type:	Rabbit IgG
Specificity:	Influenza A H7N9 (A/Shanghai/1/2013) Hemagglutinin / HA
Formulation:	0.2 µm filtered solution in PBS
Storage:	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

APPLICATIONS

Applications:	ELISA,IHC-P,FCM,ICC/IF,IF,IP
	(Antibody's applications have not been validated with corresponding viruses. Optimal concentrations/dilutions should be determined by the end user.)

RECOMMENDED CONCENTRATION

ELISA	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Influenza A H7N9 (A/Shanghai/1/2013) Hemagglutinin / HA.
--------------	---

Please Note: Optimal concentrations/dilutions should be determined by the end user.