

# 4-1BBL / CD137L / TNFSF9 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 50067-RP02

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant Mouse TNFSF9 protein (Catalog#50067-M07H)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Mouse TNFSF9 / 4-1BBL extracellular domain (rM TNFSF9; Catalog#50067-M07H; Arg 104-Glu 309; NP_033430.1). TNFSF9 specific IgG was purified by mouse TNFSF9 affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Mouse TNFSF9 / 4-1BBL
<b>Formulation:</b>	0.2 µm filtered solution in PBS with 5% trehalose
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

**Applications:** ELISA

## RECOMMENDED CONCENTRATION

**ELISA** ELISA: 1:1000-1:2000  
This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Mouse TNFSF9.

*Please Note: Optimal concentrations/dilutions should be determined by the end user.*