

# HER2 / ErbB2 Antibody, Rabbit PAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 50714-RP01

GENERAL INFORMATION	
<b>Immunogen:</b>	Recombinant Mouse HER2 / ErbB2 protein (Catalog#50714-M08H)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant Mouse HER2 / ErbB2 (rM HER2 / ErbB2; Catalog#50714-M08H; P70424; Met1-Thr653). Total IgG was purified by Protein A affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	Mouse HER2 / ErbB2
<b>Formulation:</b>	0.2 µm filtered solution in PBS with 5% trehalose
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.
APPLICATIONS	
<b>Applications:</b>	ELISA
RECOMMENDED CONCENTRATION	
<b>ELISA</b>	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect Mouse HER2 / ErbB2.

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**