

# HER2/ERBB2 Antibody, Rabbit PAb, Antigen Affinity Purified



Sino Biological  
Biological Solution Specialist

Catalog Number: 50714-T16

## GENERAL INFORMATION

|                     |   |
|---------------------|---|
| <b>Immunogen:</b>   | Recombinant Mouse HER2/ERBB2 Protein (Catalog#50714-M08H)   |
| <b>Preparation</b>  | Produced in rabbits immunized with purified, recombinant Mouse HER2/ERBB2 (rM HER2/ERBB2; Catalog#50714-M08H; P70424; Met1-Thr653). HER2/ERBB2 specific IgG was purified by Mouse HER2/ERBB2 affinity chromatography.   |
| <b>Ig Type:</b>     | Rabbit IgG  |
| <b>Specificity:</b> | Mouse HER2/ERBB2  |
| <b>Formulation:</b> | 0.2 µm filtered solution in PBS   |
| <b>Storage:</b>     | This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free.<br>Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles. |

## APPLICATIONS

**Applications:** ELISA

## RECOMMENDED CONCENTRATION

**ELISA** ELISA: 1:5000-1:10000  
This antibody can be used at 1:5000-1:10000 with the appropriate secondary reagents to detect Mouse HER2/ERBB2.

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**