

# Influenza A H5N3 Hemagglutinin / HA Antibody, Rabbit PAb



Sino Biological  
Biological Solution Specialist

Catalog Number: 11696-RP01

## GENERAL INFORMATION

<b>Immunogen:</b>	Recombinant H5N3 (A/duck/Hokkaido/167/2007) HA / Hemagglutinin protein (Catalog#11696-V08H)
<b>Preparation</b>	Produced in rabbits immunized with purified, recombinant H5N3 (A/duck/Hokkaido/167/2007) HA / Hemagglutinin (Catalog#11696-V08H; BAG07130.2; Met 1-Gln 527). Total IgG was purified by Protein A affinity chromatography.
<b>Ig Type:</b>	Rabbit IgG
<b>Specificity:</b>	H5N3 (A/duck/Hokkaido/167/2007) HA / Hemagglutinin
<b>Formulation:</b>	0.2 µm filtered solution in PBS
<b>Storage:</b>	This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C. Preservative-Free. Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## APPLICATIONS

<b>Applications:</b>	WB,ELISA
	IHC, FCM, IF, IP et al. applications haven't been validated. (Antibody's applications haven't been validated with corresponding virus positive samples. Optimal concentrations/dilutions should be determined by the end user.)

## RECOMMENDED CONCENTRATION

<b>Western Blot</b>	This antibody can be used at 1:500-1:1000 with the appropriate secondary reagents to detect H5N3 (A/duck/Hokkaido/167/2007) HA in WB.
<b>ELISA</b>	ELISA: 1:1000-1:2000 This antibody can be used at 1:1000-1:2000 with the appropriate secondary reagents to detect H5N3 (A/duck/Hokkaido/167/2007) HA.

**Please Note: Optimal concentrations/dilutions should be determined by the end user.**